



# ResMed

Actual Size:  
148 mm (Width)  
x 105 mm (Height)



Ruler is for scaling purposes.


How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

## AirFit™ N30i

Nasal cradle mask







**ResMed**

**AirFit™ N30i**  
Nasal mask

**AirFit™ P30i**  
Nasal pillows mask

## Frame Fitting Instructions



 Small frame	 Standard frame
<b>0-150mm</b>	<b>&gt;150mm</b>
<b>0-6"</b>	<b>&gt;6"</b>